



NEW MEXICO'S HEALTH INFORMATION EXCHANGE

This reference is provided at the request of Labcorp to clarify the difference in terms used in the SYNCRONYS HIE clinical portal versus the terms used by Labcorp for abnormal results.

1

SYNCRONYS/LABCORP CROSSWALK OF ABNORMAL RESULT TERMS



SYNCRONYS		Labcorp	
Abnormal Flag Value	Defined by HL7 Standard version 2.5.1 - Displayed Verbiage	Abnormal Flag Value	Defined by Labcorp External EDI Spec version 10.15
A	Abnormal	A	Abnormal
L	Below low normal	L	Below low normal
H	Above high normal	H	Above high normal
LL	Below lower panic limits	LL	Alert Low
HH	Above upper panic limits	HH	Alert High
<	Below absolute low-off instrument scale	<	Panic Low
>	Above absolute high-off instrument scale	>	Panic High
R	Resistant	R	Resistant
I	Intermediate	I	Intermediate
S	Susceptible	S	Susceptible
NEG	Negative	NEG	Negative
POS	Positive	POS	Positive

2